




## Cell collection





### Materials Required

Cold PBS	
2 pre-chilled 50ml Falcon tubes	
2 re-usable 50ml Falcon tubes labeled BALANCE	

1. Label your 2 sample 50ml Falcon tubes (white top) with the sample IDs
2. Add \_\_\_\_\_mLs (from above) of your cell culture to a pre-chilled 50-ml (**white** top) Falcon tube
3. Add \_\_\_\_\_mLs (from above) of water to a **50-ml Falcon Tube (yellow)** labeled BALANCE
4. Centrifuge at 3,000g for 5 minutes at 4°C.
5. Remove as much supernatant as possible and dispose of it into a container with 10% bleach
6. Resuspend the pellet in 10 mL ice-cold **PBS (Blue 50ml Falcon Tube)**
7. Centrifuge at 3,000g for 5 minutes at 4°C.
8. Remove as much supernatant as possible and dispose of it into a container with 10% bleach









### STEP 3 - Spheroplast (Cells without a cell wall) collection

#### Materials Required

Cold Spheroplast Buffer: 1M sorbitol, 1X PBS pH 7.4, 0.1 M EDTA	
Zymolyase (5U/ $\mu$ l)	
Cold PBS	
30°C incubator	
1.5 ml microcentrifuge tube	

1. Label your 1.5 ml tubes with your sample IDs
2. Add 500  $\mu$ L of **Spheroplast Buffer (Red 15ml Falcon Tube)** to your yeast pellet
3. Gently resuspend the pellet
4. Transfer the resuspended pellet to 1.5ml microcentrifuge tube
5. Add **10  $\mu$ l Zymolyase (5U/ $\mu$ l)] (White dot 1.5ml tube)**
6. Pipette up and down to mix with a wide-bore tip to mix.
7. Incubate at 30°C for 20 minutes.
8. Invert and gently flick tubes every 5 minutes to keep cells suspended in solution.
9. Gently pellet spheroplasts by centrifugation for 3 minutes at 300 x g
10. Remove the supernatant.
11. Using a 1000  $\mu$ l wide-bore pipette tip, gently resuspend the spheroplast pellet in 200  $\mu$ l of cold **PBS (Blue 50ml Falcon Tube)**

STEP 4 - Quick-DNA™ Miniprep Plus Kit

Resuspended Proteinase K	
BioFluid & Cell Buffer (Red)	
Genomic Binding Buffer	
DNA Pre-Wash Buffer	
g-DNA Wash Buffer	
DNA Elution Buffer	
1x Zymo-Spin column per sample	
2 Collection Tubes column per sample	
1 microcentrifuge column per sample	
55°C incubator	



1. Add 200  $\mu$ l **BioFluid & Cell Buffer (Red) (red tube)** to the sample
2. Add 20  $\mu$ l **Proteinase K (orange dot tube)** to the sample
3. Vortex for 10-15 seconds
4. Incubate at 55°C for 10 minutes
5. Add 400  $\mu$ l of **Genomic Binding Buffer (blue dot tube)**
  
6. Place a Zymo-Spin Column in a collection tube
7. Transfer 800  $\mu$ l of your mixture to the column
8. Centrifuge at 12,000 x g for 1 minute
9. Discard the flow-through and collection tube
10. Place spin column in a new collection tube
  
11. Add 400  $\mu$ l of **DNA Pre-Wash Buffer (green dot tube)** to the spin column
12. Centrifuge at 12,000 x g for 1 minute
13. Empty the collection tube
  
14. Add 700  $\mu$ l **g-DNA Wash Buffer (orange cap 15ml tube)** to the spin column.
15. Centrifuge at 12,000 x g for 1 minute
16. Empty the collection tube
  
17. Add 200  $\mu$ l **g-DNA Wash Buffer (orange cap 15ml tube)** to the spin column
18. Centrifuge at 12,000 x g for 1 minute
19. Discard the flow-through and collection tube
  
20. Transfer the spin column to a clean microcentrifuge tube
21. Add 50  $\mu$ l **DNA Elution Buffer (yellow tube)**
22. Incubate for 5 minutes at room temperature
23. Centrifuge at 12,000 x g for 1 minute
24. Your DNA will now be at the bottom of your microcentrifuge tube